


Product datasheet for **GA209619**

Mouse Poldip3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Poldip3 CRISPRa kit - CRISPR gene activation of mouse polymerase (DNA-directed), delta interacting protein 3
Symbol:	Poldip3
UniProt ID:	Q8BG81
Synonyms:	1110008P04Rik; AA408269; AL022852; C77954; mKIAA1649; PDIP46; SKAR
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)