


Product datasheet for **GA209625**

Mouse Ankrd45 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Ankrd45 CRISPRa kit - CRISPR gene activation of mouse ankyrin repeat domain 45
Symbol:	Ankrd45
Components:	GA209625G1, Ankrd45 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA209625G2, Ankrd45 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA209625G3, Ankrd45 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q810N6</u>
Synonyms:	4933409K03Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)