

Product datasheet for **GA209628**

Mouse D3ErtD751e activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	D3ErtD751e CRISPRa kit – CRISPR gene activation of mouse DNA segment, Chr 3, ERATO Doi 751, expressed
Symbol:	D3ErtD751e
Components:	GA209628G1, D3ErtD751e gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GEI00074) GA209628G2, D3ErtD751e gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GEI00074) GA209628G3, D3ErtD751e gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GEI00074) 1 CRISPRa-Enhancer vector, SKU GEI00056 1 CRISPRa scramble vector, SKU GEI00077
UniProt ID:	<u>Q8BGN2</u>
Synonyms:	2810009O15Rik; 4930415G15Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

