

Product datasheet for **GA209645**

Mouse *Irak3* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	<i>Irak3</i> CRISPRa kit - CRISPR gene activation of mouse interleukin-1 receptor-associated kinase 3
Symbol:	<i>Irak3</i>
UniProt ID:	Q8K4B2
Synonyms:	4833428C18Rik; AI563835; IRAK-M
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)