


Product datasheet for **GA209646**

Mouse Lymr1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Lymr1 CRISPRa kit - CRISPR gene activation of mouse LYR motif containing 1
Symbol:	Lymr1
Components:	GA209646G1, Lymr1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA209646G2, Lymr1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA209646G3, Lymr1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9CQB7</u>
Synonyms:	1110065L10Rik; 2310004B22Rik; 4930404J24Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)