

Product datasheet for **GA209652**

Mouse Hapln2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Hapln2 CRISPRa kit - CRISPR gene activation of mouse hyaluronan and proteoglycan link protein 2
Symbol:	Hapln2
Components:	GA209652G1, Hapln2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA209652G2, Hapln2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA209652G3, Hapln2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9ESM3</u>
Synonyms:	4930401E20Rik; Bral1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)