

Product datasheet for **GA209716**

Mouse Slc35a5 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Slc35a5 CRISPRa kit - CRISPR gene activation of mouse solute carrier family 35, member A5
Symbol:	Slc35a5
Components:	GA209716G1, Slc35a5 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA209716G2, Slc35a5 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA209716G3, Slc35a5 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q921R7</u>
Synonyms:	1010001J06Rik; AU021179; BB097433; D16Ert450e; D730043G07Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)