

Product datasheet for **GA209724**

Mouse Pi16 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Pi16 CRISPRa kit - CRISPR gene activation of mouse peptidase inhibitor 16
Symbol:	Pi16
Components:	GA209724G1, Pi16 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA209724G2, Pi16 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA209724G3, Pi16 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9ET66</u>
Synonyms:	1200009H11Rik; Cripi; PI-16
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)