

Product datasheet for **GA209734**

Mouse Sash3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Sash3 CRISPRa kit - CRISPR gene activation of mouse SAM and SH3 domain containing 3
Symbol:	Sash3
Components:	GA209734G1, Sash3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA209734G2, Sash3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA209734G3, Sash3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8K352</u>
Synonyms:	1200013B08Rik; AW413946; SLY1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)