


Product datasheet for **GA209757**

Mouse Tmem38a activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tmem38a CRISPRa kit - CRISPR gene activation of mouse transmembrane protein 38A
Symbol:	Tmem38a
Components:	GA209757G1, Tmem38a gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA209757G2, Tmem38a gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA209757G3, Tmem38a gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q3TMP8</u>
Synonyms:	1110001E17Rik; AI413399; mg33a; TRIC-A
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)