

Product datasheet for **GA209774**

Mouse P2ry13 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	P2ry13 CRISPRa kit - CRISPR gene activation of mouse purinergic receptor P2Y, G-protein coupled 13
Symbol:	P2ry13
UniProt ID:	<u>Q9D8I2</u>
Synonyms:	2010001L06Rik; GPCR1; Gpr86; GPR94; P2Y13; P2yr13; SP174
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)