

Product datasheet for **GA209789**

Mouse *Tex37* activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | Tex37 CRISPRa kit - CRISPR gene activation of mouse testis expressed 37 |
| Symbol: | Tex37 |
| UniProt ID: | Q9DAG4 |
| Synonyms: | 1700011F03Rik |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)