

Product datasheet for **GA209798**

Mouse Slx4ip activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Slx4ip CRISPRa kit - CRISPR gene activation of mouse SLX4 interacting protein
Symbol:	Slx4ip
Components:	GA209798G1, Slx4ip gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA209798G2, Slx4ip gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA209798G3, Slx4ip gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9D7Y9</u>
Synonyms:	2210009G21Rik; 2410004I22Rik; C85661; Slxip
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)