

Product datasheet for **GA209830**

Mouse Wdr33 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Wdr33 CRISPRa kit - CRISPR gene activation of mouse WD repeat domain 33
Symbol:	Wdr33
UniProt ID:	<u>Q8K4P0</u>
Synonyms:	1110001N06Rik; 2310011G05Rik; 2810021O11Rik; 8430413N20Rik; WDC146
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)