


Product datasheet for **GA209852**

Mouse Lonrf3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Lonrf3 CRISPRa kit - CRISPR gene activation of mouse LON peptidase N-terminal domain and ring finger 3
Symbol:	Lonrf3
UniProt ID:	Q9D4H7
Synonyms:	4932412G04Rik; 5730439E01Rik; A830039N02Rik; AU023707; Rnf127
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)