


Product datasheet for **GA209884**

Mouse **Lrriq3** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	Lrriq3 CRISPRa kit - CRISPR gene activation of mouse leucine-rich repeats and IQ motif containing 3
Symbol:	Lrriq3
Components:	GA209884G1, Lrriq3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA209884G2, Lrriq3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA209884G3, Lrriq3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q14DL3</u>
Synonyms:	4930438B07Rik; 4930511J15Rik; 4933403H06Rik; BB020367; Lrrc44
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)