


## Product datasheet for **GA209891**

### Mouse Arl13a activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                               |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                                  |
| Description:  | Arl13a CRISPRa kit - CRISPR gene activation of mouse ADP-ribosylation factor-like 13A                                                                                                                                                                                                                         |
| Symbol:       | Arl13a                                                                                                                                                                                                                                                                                                        |
| Components:   | GA209891G1, Arl13a gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA209891G2, Arl13a gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA209891G3, Arl13a gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <u><a href="#">Q9D416</a></u>                                                                                                                                                                                                                                                                                 |
| Synonyms:     | 4933422M21Rik; Arl13                                                                                                                                                                                                                                                                                          |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                       |

#### Product images:



[View online »](#)