

## Product datasheet for **GA209907**

### Mouse Snx29 activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                            |
|---------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                               |
| Description:  | Snx29 CRISPRa kit - CRISPR gene activation of mouse sorting nexin 29                                                                                                                                                                                                                                       |
| Symbol:       | Snx29                                                                                                                                                                                                                                                                                                      |
| Components:   | GA209907G1, Snx29 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA209907G2, Snx29 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA209907G3, Snx29 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <u><a href="#">Q9D3S3</a></u>                                                                                                                                                                                                                                                                              |
| Synonyms:     | 4933437K13Rik; ENSMUSG00000079738; Gm930; Gm1737; Gm11170; Rundc2a                                                                                                                                                                                                                                         |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                    |

#### Product images:



[View online »](#)