

Product datasheet for **GA209928**

Mouse Fam234b activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Fam234b CRISPRa kit - CRISPR gene activation of mouse family with sequence similarity 234, member B
Symbol:	Fam234b
Components:	GA209928G1, Fam234b gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA209928G2, Fam234b gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA209928G3, Fam234b gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8BYI8</u>
Synonyms:	8430419L09Rik; Kiaa1467; mKIAA1467
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)