

Product datasheet for **GA209937**

Mouse Nipal3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Nipal3 CRISPRa kit - CRISPR gene activation of mouse NIPA-like domain containing 3
Symbol:	Nipal3
Components:	GA209937G1, Nipal3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA209937G2, Nipal3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA209937G3, Nipal3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8BGN5</u>
Synonyms:	9130020G22Rik; AI480660; Npal3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)