

Product datasheet for **GA209960**

Mouse *Scrn3* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Scrn3 CRISPRa kit - CRISPR gene activation of mouse secernin 3
Symbol:	Scrn3
Components:	GA209960G1, Scrn3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA209960G2, Scrn3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA209960G3, Scrn3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q3TMH2</u>
Synonyms:	4833415E20Rik; AI427019; AV002033; Ses3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)