

## Product datasheet for **GA209974**

### Mouse Brme1 activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                            |
|---------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                               |
| Description:  | 4930432K21Rik CRISPRa kit - CRISPR gene activation of mouse RIKEN cDNA 4930432K21 gene                                                                                                                                                                                                                     |
| Symbol:       | Brme1                                                                                                                                                                                                                                                                                                      |
| Components:   | GA209974G1, Brme1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA209974G2, Brme1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA209974G3, Brme1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <u><a href="#">Q6DIA7</a></u>                                                                                                                                                                                                                                                                              |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                    |

#### Product images:



[View online »](#)