

Product datasheet for **GA209978**

Mouse Ptchd3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Ptchd3 CRISPRa kit - CRISPR gene activation of mouse patched domain containing 3
Symbol:	Ptchd3
Components:	GA209978G1, Ptchd3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA209978G2, Ptchd3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA209978G3, Ptchd3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q0EEE2</u>
Synonyms:	4930451E13Rik; 4933440L20Rik; RNDEu-3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)