

Product datasheet for **GA209979**

Mouse Wdr35 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Wdr35 CRISPRa kit - CRISPR gene activation of mouse WD repeat domain 35
Symbol:	Wdr35
Components:	GA209979G1, Wdr35 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA209979G2, Wdr35 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA209979G3, Wdr35 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8BND3</u>
Synonyms:	4930459M12Rik; 4931430C06; mKIAA1336
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)