

Product datasheet for **GA209981**

Mouse *Lrrc74b* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Lrrc74b CRISPRa kit – CRISPR gene activation of mouse leucine rich repeat containing 74B
Symbol:	Lrrc74b
Components:	GA209981G1, Lrrc74b gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA209981G2, Lrrc74b gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA209981G3, Lrrc74b gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q14BP6</u>
Synonyms:	4930451C15Rik; mCG12354; mCP17390
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

