


## Product datasheet for **GA210003**

### Mouse C2cd5 activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                            |
|---------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                               |
| Description:  | C2cd5 CRISPRa kit - CRISPR gene activation of mouse C2 calcium-dependent domain containing 5                                                                                                                                                                                                               |
| Symbol:       | C2cd5                                                                                                                                                                                                                                                                                                      |
| Components:   | GA210003G1, C2cd5 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA210003G2, C2cd5 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA210003G3, C2cd5 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <u><a href="#">Q7TPS5</a></u>                                                                                                                                                                                                                                                                              |
| Synonyms:     | 5730419I09Rik; C030008B15Rik; CDP138                                                                                                                                                                                                                                                                       |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                    |

#### Product images:



[View online »](#)