

Product datasheet for **GA210013**

Mouse **Pik3cb** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	Pik3cb CRISPRa kit - CRISPR gene activation of mouse phosphatidylinositol-4,5-bisphosphate 3-kinase catalytic subunit beta
Symbol:	Pik3cb
UniProt ID:	<u>Q8BTI9</u>
Synonyms:	1110001J02Rik; AI447572; p110beta
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)