

Product datasheet for **GA210034**

Mouse Morn3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Morn3 CRISPRa kit - CRISPR gene activation of mouse MORN repeat containing 3
Symbol:	Morn3
Components:	GA210034G1, Morn3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA210034G2, Morn3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA210034G3, Morn3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8C5T4</u>
Synonyms:	4930438O03Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)