

## Product datasheet for GA210068

### Mouse Tmbim7 activation kit by CRISPRa

#### Product data:

|                      |                                                                                                                                                                                                                                                                                                               |
|----------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type:</b> | CRISPRa Kits                                                                                                                                                                                                                                                                                                  |
| <b>Description:</b>  | Tmbim7 CRISPRa kit – CRISPR gene activation of mouse transmembrane BAX inhibitor motif containing 7                                                                                                                                                                                                           |
| <b>Symbol:</b>       | Tmbim7                                                                                                                                                                                                                                                                                                        |
| <b>Components:</b>   | GA210068G1, Tmbim7 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA210068G2, Tmbim7 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA210068G3, Tmbim7 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| <b>Synonyms:</b>     | 4930403J02Rik; 4930500J03Rik; 4930511M11Rik; Lfg5; Lgf5; Tmbim1b                                                                                                                                                                                                                                              |
| <b>Format:</b>       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                       |

#### Product images:

