

Product datasheet for **GA210077**

Mouse Purg activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Purg CRISPRa kit - CRISPR gene activation of mouse purine-rich element binding protein G
Symbol:	Purg
Components:	GA210077G1, Purg gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA210077G2, Purg gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA210077G3, Purg gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8R4E6</u>
Synonyms:	4930486B15Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)