

Product datasheet for **GA210095**

Mouse Uhrf1bp1l activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Uhrf1bp1l CRISPRa kit - CRISPR gene activation of mouse UHRF1 (ICBP90) binding protein 1-like
Symbol:	Uhrf1bp1l
Components:	GA210095G1, Uhrf1bp1l gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA210095G2, Uhrf1bp1l gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA210095G3, Uhrf1bp1l gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	A2RSJ4
Synonyms:	2010319N22Rik; 4930506D01Rik; AA589495; C80672; E030041M21Rik; mKIAA0701
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)