

Product datasheet for **GA210122**

Mouse Tex38 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tex38 CRISPRa kit - CRISPR gene activation of mouse testis expressed 38
Symbol:	Tex38
Components:	GA210122G1, Tex38 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA210122G2, Tex38 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA210122G3, Tex38 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>A2A8T7</u>
Synonyms:	4930544O15Rik; Theg4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)