

Product datasheet for **GA210131**

Mouse Prr3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Prr3 CRISPRa kit - CRISPR gene activation of mouse proline-rich polypeptide 3
Symbol:	Prr3
UniProt ID:	Q811B5
Synonyms:	4930540G07Rik; Cat56
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)