

Product datasheet for **GA210159**

Mouse Prkd3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Prkd3 CRISPRa kit - CRISPR gene activation of mouse protein kinase D3
Symbol:	Prkd3
UniProt ID:	Q8K1Y2
Synonyms:	4930557O20Rik; 5730497N19Rik; Pkcnu; PKD3; Prkcn
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)