

Product datasheet for GA210169

Mouse Taf1d activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Taf1d CRISPRa kit - CRISPR gene activation of mouse TATA-box binding protein associated factor, RNA polymerase I, D
Symbol:	Taf1d
Components:	GA210169G1, Taf1d gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA210169G2, Taf1d gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA210169G3, Taf1d gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9D4V4</u>
Synonyms:	2810003M17Rik; 4930553M18Rik; 5930426I11Rik; Josd3; TAF(I)41; TAFI41
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

