

Product datasheet for **GA210175**

Mouse Fam243 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Fam243 CRISPRa kit – CRISPR gene activation of mouse family with sequence similarity 243
Symbol:	Fam243
Components:	GA210175G1, Fam243 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA210175G2, Fam243 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA210175G3, Fam243 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8CDS7</u>
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

