

Product datasheet for **GA210199**

Mouse Arhgap36 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Arhgap36 CRISPRa kit - CRISPR gene activation of mouse Rho GTPase activating protein 36
Symbol:	Arhgap36
Components:	GA210199G1, Arhgap36 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA210199G2, Arhgap36 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA210199G3, Arhgap36 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>B1AUC7</u>
Synonyms:	1100001E04Rik; rho-gap
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)