


Product datasheet for **GA210247**

Mouse Tex48 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tex48 CRISPRa kit - CRISPR gene activation of mouse testis expressed 48
Symbol:	Tex48
Components:	GA210247G1, Tex48 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA210247G2, Tex48 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA210247G3, Tex48 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	bM230C2.3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)