

Product datasheet for **GA210257**

Mouse Paqr9 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Paqr9 CRISPRa kit - CRISPR gene activation of mouse progesterin and adipoQ receptor family member IX
Symbol:	Paqr9
Components:	GA210257G1, Paqr9 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA210257G2, Paqr9 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA210257G3, Paqr9 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6TCG2</u>
Synonyms:	1700020G04Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)