

Product datasheet for GA210262

Mouse Ep400 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Ep400 CRISPRa kit - CRISPR gene activation of mouse EIA binding protein p400
Symbol:	Ep400
Components:	GA210262G1, Ep400 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA210262G2, Ep400 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA210262G3, Ep400 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8CHI8</u>
Synonyms:	1700020J09Rik; AU023439; mDomino; mKIAA1498; p400
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

