

Product datasheet for **GA210290**

Mouse Tex30 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tex30 CRISPRa kit - CRISPR gene activation of mouse testis expressed 30
Symbol:	Tex30
Components:	GA210290G1, Tex30 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA210290G2, Tex30 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA210290G3, Tex30 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q3TUU5</u>
Synonyms:	1700029F09Rik; 3110030D08Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)