


Product datasheet for **GA210302**

Mouse Lin37 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Lin37 CRISPRa kit - CRISPR gene activation of mouse lin-37 homolog (C. elegans)
Symbol:	Lin37
Components:	GA210302G1, Lin37 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA210302G2, Lin37 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA210302G3, Lin37 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9D8N6</u>
Synonyms:	1810054G18Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)