

Product datasheet for **GA210304**

Mouse Tmem225 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tmem225 CRISPRa kit - CRISPR gene activation of mouse transmembrane protein 225
Symbol:	Tmem225
Components:	GA210304G1, Tmem225 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA210304G2, Tmem225 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA210304G3, Tmem225 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9D9S2</u>
Synonyms:	1700030E15Rik; Tmem225
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)