

Product datasheet for GA210306

Mouse Tex22 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tex22 CRISPRa kit - CRISPR gene activation of mouse testis expressed gene 22
Symbol:	Tex22
Components:	GA210306G1, Tex22 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA210306G2, Tex22 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA210306G3, Tex22 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9D9U4
Synonyms:	1700028O09Rik; AV210422; Tep22
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)