

Product datasheet for **GA210333**

Mouse Svip activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Svip CRISPRa kit - CRISPR gene activation of mouse small VCP/p97-interacting protein
Symbol:	Svip
Components:	GA210333G1, Svip gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA210333G2, Svip gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA210333G3, Svip gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q3UZP4</u>
Synonyms:	1700006C13Rik; 6620401M08Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)