

Product datasheet for **GA210343**

Mouse Plppr5 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Plppr5 CRISPRa kit - CRISPR gene activation of mouse phospholipid phosphatase related 5
Symbol:	Plppr5
Components:	GA210343G1, Plppr5 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA210343G2, Plppr5 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA210343G3, Plppr5 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q8BJ52
Synonyms:	4833424O15Rik; Lppr5; Pap2d; PRG-5
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)