


Product datasheet for **GA210347**

Mouse Ttc23l activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Ttc23l CRISPRa kit - CRISPR gene activation of mouse tetratricopeptide repeat domain 23-like
Symbol:	Ttc23l
Components:	GA210347G1, Ttc23l gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA210347G2, Ttc23l gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA210347G3, Ttc23l gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	A6H6E9
Synonyms:	4930401A09Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)