

Product datasheet for **GA210359**

Mouse Tasp1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tasp1 CRISPRa kit - CRISPR gene activation of mouse taspase, threonine aspartase 1
Symbol:	Tasp1
UniProt ID:	<u>Q8R1G1</u>
Synonyms:	4930485D02Rik; AW986064
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)