

Product datasheet for **GA210412**

Mouse Dock11 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Dock11 CRISPRa kit - CRISPR gene activation of mouse dedicator of cytokinesis 11
Symbol:	Dock11
Components:	GA210412G1, Dock11 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA210412G2, Dock11 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA210412G3, Dock11 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>A2AF47</u>
Synonyms:	5033414A21Rik; 8030476J24; sph238; Zizimin2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)