


Product datasheet for **GA210419**

Mouse **Gon4l** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	Gon4l CRISPRa kit - CRISPR gene activation of mouse gon-4-like (C.elegans)
Symbol:	Gon4l
Components:	GA210419G1, Gon4l gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA210419G2, Gon4l gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA210419G3, Gon4l gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	1500041I23Rik; 2610100B20Rik; 3632445O20Rik; 5430404N14Rik; AA387005
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)